


BX2WI SYSTEM DIAGRAM



*1 Slight vignetting may occur in the periphery of the field of view in combination with an additional intermediate attachment. *2 Slight vignetting may occur in the periphery of the field of view in combination with fluorescence illuminator. *3 Can be used with U-ETR and U-TR30IR. Field of view is limited when using more than two intermediate tubes. *4 Sub port can accept U-TV1x-2, and 2/3" or less CCD cameras. *5 F.N. is 22 when fluorescence observation. *6 Use with fluorescent illuminator. *7 Acceptable video adapters are U-TV1x-2, U-TVACA, U-PMTVC2xIR and U-PMTVC4x.